


SATURN 8000 1717N

Wide and Slim Portable Flat Panel Detector for Digital Radiography



anytime™

Stable and reliable AED
(Automatic Exposure
Detection)



inside AP™

Fast wireless
communication



High image quality




Easy installation
(cassette tray)

Features

- Wide active area of 17" x 17"
- High spatial resolution with 140um pixel array
- Stable and reliable automatic exposure detection (Anytime™)
- Fast wireless communication by IEEE 802.11n (Inside AP™)
- Viewer software running on Windows™ OS (VXvue™)
- Various applications such as General Radiography, Orthopedic, Veterinary, Chiropractic X-ray systems

Drawing



Technical Specifications

Model Name	Saturn 8000 17x17N Csl / Saturn 8000 17x17N Gadox	Saturn 8000 17x17NW Csl / Saturn 8000 17x17NW Gadox
Technology	a-Si TFT with photodiode	
Scintillator	Csl:TI	Gd ₂ O ₂ S:Tb
Pixel Pitch	140 μm	
Spatial Resolution	3.5 lp/mm	
Pixels	3072 x 3072 pixels	
Image Size	17" x 17" (43.0 cm x 43.0 cm)	
Grayscale	16 bit	
Image Acquisition Time	Wired: 1.5 s Wireless: 3 s	
Recommended Cycle Time	15 s	
Data Interface	Gigabit Ethernet IEEE 802.11n	
X-ray Generator Interface	Line trigger: DR trigger Mode Auto trigger: AED Mode	
Dimensions	46.0 cm x 46.0 cm x 1.5 cm	
Weight	4.2 kg	4.5 kg (including battery)
Battery	-	Lithium Ion 3100 mAh x 2 6 h (Capturing)
Dust and Water Resistant	IP 53	
X-ray Voltage Range	40 – 150 kVp	
Operating Environment	10°C to 35°C 30 % to 85 % RH (non-condensing)	
Power	DC 24V, 1.0A	

* Specifications are subject to change without prior notice.

Office

3F., NO.32, Sec. 2, Zhongyong S. Rd., Beitou Dist., Taipei City 11270, Taiwan
Tel +886-2-2892-2275 Fax +886-2-2892-2172 E-mail Sales@newmedical.com.tw

Factory

41-3, Burim-ro 170beon-gil, Dongan-gu, Anyang-si, Gyeonggi-do, 14055
Republic of KOREA

www.nmi-dr.com

FDA Approved  